

SPEAKING WITH ANNUAL **REFILL** PATIENTS?

Give them an option that fits in their pocket

AUVI-Q WAS DESIGNED TO BE EASY TO CARRY AND USE.



POCKET SIZED



VOICE INSTRUCTIONS



AUTO-RETRACTABLE NEEDLE

Patients should seek emergency medical care immediately after use.

AUVI-Q is available in 3 strengths:



AUVI-Q 0.1 mg
For infants and toddlers¹
weighing 16.5 to 33 lbs
(7.5 to 15 kg)



AUVI-Q 0.15 mg
For children¹
weighing 33 to 66 lbs
(15 to 30 kg)



AUVI-Q 0.3 mg
For anyone¹
weighing 66 lbs
(30 kg) or more

Train patients by phone or telehealth by using an AUVI-Q Trainer or by sending them to auvi-q.com/training to watch a step-by-step video.

LOW OUT-OF-POCKET COSTS FOR ELIGIBLE COMMERCIALLY INSURED PATIENTS

with no prior authorization required. Check patient
eligibility at auvi-q.com/savings



Scan for information about
the AUVI-Q Savings Offer

PRESCRIBE THROUGH:



IN-STORE PICKUP:

Deselect “Use generic”
Select Walgreens in EMR

Pharmacy team members can review
prescription processing details at:

StoreNet > Rx > Patient Care
and search for “Kaleo”



HOME DELIVERY:

Select ASPN

(mail order pharmacy) in EMR

Patients can text **AUVI-Q** to 84790
to manage their prescription via text*,
or they can call 1-844-357-3968 to
schedule their delivery

*Message and data rates may apply.

OR

Remember to deselect any preset preferred pharmacy for your patient's AUVI-Q prescription.

Each carton contains
2 AUVI-Q and 1 Trainer



Reference: 1. AUVI-Q. Prescribing information. kaleo, Inc.; 2019.

AUVI-Q[®] and AUVI-q[®] are registered trademarks of kaleo, Inc. Other referenced trademarks are owned by their respective third party owners.
© 2022 kaleo, Inc. CM-US-AQ-2582